

Quality grades for every analytical application

Our portfolio of brands and product grades offer a range of solutions for your analytical testing needs.

- Application-specific testing for high resolution chromatography including, LC-MS, UHPLC-MS and GC Headspace
- Rigorous quality assurance and testing procedures throughout the production process, ensure the lot-to-lot consistency required for interference-free reproducible results every time
- High-volume solvent delivery systems, available in 10L to 1000L, offer environmentally friendly solvent handling solutions for your applications, enhancing safety and improving productivity within your lab

Fields	Application	Grade(s)	Description and Use
Chromatography	Liquid Chromatography (LC)	HPLC	Suitable for HPLC gradient
		Optima LC-MS	Suitable for LC-MS and UHPLC-UV
		Thermo Scientific UHPLC-MS	Highest purity for UHPLC-MS
	Gas Chromatography (GC)	Distol	For pesticides and residue analysis
		GC Headspace	Tested to ensure absence of Type 1, 2 and 3 volatiles
LC, GC, Spectrophometry & ACS	Thermo Scientific Chromplete	Universal choice for multiple techniques	
Elemental Analysis	AAS	For analysis Primar Plus	Acids analysis for metals in the very low parts per million (ppm) range
	ICP-OES	Primar Plus TraceMetal	Acids analyzed for 55 metals in low to very low parts per billion (ppb) range
	ICP-MS	Optima	Acids analyzed for more than 65 metals in the parts per trillion (ppt) range
General Analytical Techniques	Karl Fischer Titration	Aqualine	Karl Fischer titration reagents
	Volumetric titration	For volumetric analysis	Acid and base solutions suitable for precise titration
	pH titration	For pH Metry	Buffer standard solution
	Other Analysis	For Analysis Extra Pure	Purity reagents suitable for general lab techniques

thermo
scientific

F **fisher**
chemical

To place an order, contact **your local Fisher Scientific office**

Or visit eu.fishersci.com/go/food-beverage

Austria: +43(0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
 Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Germany: +49 (0)2304 9325
 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478 Lithuania: +370 5 2444442
 Netherlands: +31 (0)20 487 70 00 Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50
 Spain: +34 902 239 303 Sweden: +46 31 352 32 00 Switzerland: +41 (0)56 618 41 11 UK: +44 (0)1509 555 500

F **fisher scientific**
part of Thermo Fisher Scientific