

Replace your old sample purification system with a new Thermo Scientific KingFisher Flex platform and **get 35% off!***

**Limited
Time
Savings**

**We will
even pick up
your old system
and dispose of
it for you.****



Fast, Flexible, Reliable

The Thermo Scientific™ KingFisher™ Flex automated purification system is designed to provide quality yields of nucleic acids, proteins and cells from a variety of sources. This versatile platform offers excellent performance and convenience, and helps you keep your cost-per-run under control.

Offer ends December 31st



Replace your old purification system and get 35% off!


KingFisher Flex purification systems provide excellent performance and flexibility for a range of applications.

- **High-throughput system can process up to 96 samples from a variety of sample materials**
- **Use with optimized KingFisher Pure nucleic acid purification kits for a complete solution**
- **Open, flexible system designed to work with any magnetic bead-based applications**
- **Paired with Thermo Scientific™ BindIt™ software for convenient instrument control, protocol development and modification**



Product Code Information

Cat. No.	Description
22-387-005	KingFisher Flex with 96 deepwell head
22-387-006	KingFisher Flex with 24 deepwell head



* Please order using these special model numbers to take advantage of package pricing on purification systems and environmentally-friendly disposal of trade-in instrument. Please include your contact information with your order. You will be contacted by Thermo Fisher Scientific and sent a Bio-Hazard Certificate for Disposal form which must be completed and returned per instructions on the form, prior to the shipment of your new purification system. Once the form is received back to Thermo Fisher Scientific, and your purification system ships, you will be contacted by the carrier to schedule the delivery date. In the event that your trade-in instrument is not ready to be removed at the time of the scheduled delivery date and a return trip by the carrier to pick up your trade-in instrument is necessary, a \$526.00 fee will apply. ***Instrument trade-in required to take advantage of this limited time promotional pricing offer. Please have trade-in instrument asset tags removed before the scheduled delivery date. **Instrument decontamination form must be received by Thermo Fisher Scientific before it can be removed from site.**

Terms and conditions: Valid in the United States and Canada for orders placed between July 1, 2015 and December 31, 2015, and only in connection with the product numbers listed above. "Savings" represents the difference between the manufacturer's suggested list price for individual purchase of each item listed above and the manufacturer's suggested promotional price. All offers subject to Thermo Fisher standard terms and conditions of sale. Subject to dealer participation. Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. Items shown at no cost represent "in-kind discounts" and the value of the specific offer(s) the customer may receive from the manufacturer under this program is a "discount or other reduction in price" to the customer under section 1128B(b)(3)(A) of the Social Security Act [42 U.S.C. 1320a-7b (b)(3)(a)]. Accordingly, the customer shall disclose this and any other discounts or other reduction in price received under this program under any state or federal program which provides cost or charge-based reimbursements to the customer for the products or services purchased under this program.

© 2015 Thermo Fisher Scientific. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries.

www.thermoscientific.com/kingfisher-offer

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries.



In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: www.fishersci.com

PFKFFLEXNA 0815 BNO819153

In Canada:

For customer service, call 1-800-234-7437

To fax an order, use 1-800-463-2996

To order online: www.fishersci.ca

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand